



SonixCAM

The easy way to record complete ultrasound procedures



SonixCAM captures entire ultrasound procedures on video for training and documentation

SonixCAM is an optional accessory for Sonix Series Systems. The adjustable HD Camera connects to the system to stream and record high-quality videos. Using a picture-in-picture feature, SonixCAM captures both the ultrasound images and video from the HD camera in one synchronized file.

Why SonixCAM?



Document your work

Self produce your own videos without expensive production equipment or additional personnel.



Produce educational videos

Record videos of your department's scanning procedures to simplify and standardize training.





Live Telemedicine

Stream live scanning to colleagues anywhere in the world over the internet.

Customize each video

SonixCAM eliminates the need for post production editing to identify your videos. Using SonixCAM software, it's easy to add title screens before and after each video.



Common format for online sharing

SonixCAM records videos in MPEG-1 movie format, which is easy to email to colleagues. Or you can upload the videos to popular sharing sites.



For more information:
Call 1-866-437-9508 or visit www.ultrasonix.com/camera